

User Manual for AutoMATA

1. Introduction	2
2. Environment	2
3. The workflow of AutoMATA	3
3.1 Model training	4
3.2 Model application	6
3.3 Data analysis	8
1: Correlation Heatmap	9
2: differential gene cluster heatmap	10
3: PCA	12
4: VENN	13
5: volcano with GSEA	14
6: GO Enrichment	15
7: KEGG Enrichment	16
8: Dumbbell_Bar	17
9: DESeq2 for read count	19
10: limma for fpkm	19

1. Introduction

AutoMATA is a user-friendly analysis platform designed for the processing and analysis of gene, mRNA, and protein expression data. Tailored for multi-omics studies in *Homo sapiens*, *Mus musculus*, *Bos taurus*, and *Drosophila melanogaster*, AutoMATA streamlines the analysis pipeline from raw data to biological insights. The platform integrates eight deep learning models alongside sixteen analytical modules, providing a flexible framework for model training, interpretation, and downstream biological analysis. AutoMATA is organised into three major components: (A) the Data Processing Module, (B) the Deep Learning Module, and (C) the Data Analysis Module. Together, these three modules form a cohesive workflow that bridges multi-omics data processing, predictive modelling, and bioinformatics analysis. AutoMATA is freely available at <http://automata.biotoools.bio/>, and its source code is provided on GitHub at <https://github.com/ABILiLab/AutoMATA>.

2. Environment

To simplify the installation process and reduce the complexity in environment configuration, AutoMATA provides pre-packaged environment files in .yaml format. Users can set up the required environments by executing a single command in the terminal.

AutoMATA's runtime environment is divided into two parts based on its functional modules:

- Model Training and Application Environment

This environment is developed based on Python and includes essential packages such as PyTorch (1.10.1), scikit-learn (0.24.2), pandas (1.1.5), and numpy (1.19.2).

To install this environment, use the following command:

```
conda env create -f environment.yaml
```

- Data Analysis Environment

This environment is based on the R language. The environment configuration file is "environment_R.yaml", and the installation command is:

```
conda env create -f environment_R.yaml
```

Then, there are also some R packages that need to be installed by running R scripts. First, activate your configured conda environment:

```
conda activate env_name
```

(Replace “env_name” with the name of your environment.)

Run the script ‘install_packages.R’:

```
Rscript install_packages.R
```

3. The workflow of AutoMATA

This section provides a step-by-step guide to the overall workflow of AutoMATA, which is organised into three models: (A) the Data Processing Module, (B) the Deep Learning Module, and (C) the Data Analysis Module.

The Deep Learning Module incorporates eight deep learning (DL) models for both model training and application. These models include:

- Convolutional Neural Network (CNN)
- Autoencoder
- Long Short-Term Memory (LSTM)
- Multilayer Perception (MLP)
- Recurrent Neural Network (RNN)
- Radial Basis Function Neural Network (RBFNN)
- Self-Organising Map (SOM)
- Transformer

The Data Analysis Module provides a range of bioinformatics analysis functions in omics data implemented in R. This module encompasses a variety of functions for differential expression analysis, data visualisation, and functional enrichment. Key features include:

- Correlation heatmap
- Differential gene cluster heatmap
- Principal component analysis (PCA)
- Venn diagrams (VENN)
- Volcano plots with Gene Set Enrichment Analysis (GSEA)
- Gene Ontology (GO) enrichment analysis
- Kyoto Encyclopedia of Genes and Genomes (KEGG) pathway enrichment
- Dumbbell and bar plots (Dumbbell_Bar)
- DESeq2 for read count-based expression data
- limma for FPKM-based expression analysis

3.1 Model training

For the Deep Learning Module in AutoMATA, we provide two flexible training strategies to accommodate different research needs and ensure robust model performance: the traditional train-validation-test split and stratified k-fold cross-validation. Users can select their preferred approach, upload the corresponding data, and customise the training process to suit their specific analytical goals. AutoMATA allows fine-tuning several key hyperparameters, including the number of epochs, learning rate, early stopping patience, and batch size. Additionally, users can specify the loss function and optimiser based on the distribution and nature of their dataset. The number of output classes can be set between 2 and 7 to support binary and multi-class classification tasks.

Environment Setup

To begin model training, users must first configure the Python-based conda environment as described in the “Environment” section. After activating the environment, navigate to the appropriate directory and execute the training script via the command line using the desired parameters.

Script Parameters

- `--kfold`: Integer; set to 0 to disable k-fold, or to a number such as 3 to enable stratified k-fold cross-validation with three folds.
- `--ratio`: String; defines dataset split (e.g., "8:1:1" for train/validation/test). Set to 0 to disable.
- `--epochs`: Number of training epochs.
- `--es`: Early stopping patience. Set equal to epochs to disable early stopping.
- `--lr`: Learning rate.
- `--bs`: Batch size.
- `--loss_function`: Choose from crossentropy, nllloss, or focalloss.
- `--optimizer_function`: Choose from adam, rmsprop, or sgd.
- `--output_size`: Number of output classes, the range is [2, 7].




Figure 1. Example of accuracy-loss curves of the RBFNN model after hyperparameter tuning.

Data Format

The training data should be in tab-separated (.txt) format. The first column must contain sample names, and the last column should be labelled "Label" and contain integer class labels ranging from 0 to 5. The first row must contain column headers.

Example Commands

Activate your configured conda environment:

```
conda activate env_name
```

(Replace “env_name” with the name of your environment.)

Navigate to the training directory:

```
cd code/train
```

Run model training using the following commands:

- Autoencoder

```
python autoencoder.py --kfold 0 --ratio 0 --epochs 50 --es 10 --lr 0.01 --bs 32  
--loss_function crossentropy --optimizer_function adam --output_size 2
```

- CNN

```
python cnn.py --kfold 0 --ratio 8:1:1 --epochs 50 --es 10 --lr 0.01 --bs 32  
--loss_function crossentropy --optimizer_function adam --output_size 4
```

- LSTM

```
python lstm.py --kfold 4 --ratio 0 --epochs 50 --es 10 --lr 0.005 --bs 32  
--loss_function nllloss --optimizer_function rmsprop --output_size 2
```

- MLP

```
python mlp.py
```

You may also adjust internal model parameters such as “hidden_size_1” to match the distribution of your dataset, or “dropout_rate” to prevent overfitting.

3.2 Model application

The Model Application functionality within the Deep Learning Module enables users to apply trained models to predict unknown samples – an essential step for downstream biological discovery. Building upon the previously trained models, users can select a model type, upload a trained model, and provide a testing dataset. AutoMATA will then generate two output files: one reporting the model’s testing performance and another containing predicted labels for unknown samples, thereby

addressing predictive needs in multi-omics research.

Usage Instructions

After completing model training and setting up the appropriate conda environment (as described in the “Model training” section), users can apply their trained models using the script commands outlined below.

Script Parameters

- --bs: Batch size.
- --model_type: Specify one of the supported model types: CNN, AutoEncoder, LSTM, MLP, RBFN, RNN, or Transformer.
- --model_path: Path to the trained model file (e.g., model.pth).
- --model_autoencoder_path: Required only when the ‘model_type’ is ‘AutoEncoder’; specify the path to the saved ‘model_autoencoder.pth’ file.

Input Data Format

The input test file should be tab-delimited (.txt) with:

- The first column: Sample names
- The last column: Ground truth labels (integer values ranging from 0 to 5)
- The first row: Column headers; the last column must be named “Label”

Example Commands

1. Activate your conda environment:

```
conda activate env_name
```

(Replace ‘env_name’ with your environment name.)

2. Navigate to the application directory:

```
cd code/application
```

3. Run the model application:

To apply CNN (default parameters):

```
python general.py
```

To apply the Transformer:

```
python general.py --model_type Transformer --model_path ./model/model.pth --bs  
16
```

To apply the AutoEncoder:

```
python general.py --model_type AutoEncoder  
--model_path ./model/model_cls.pth  
--model_autoencoder_path ./model/model_autoencoder.pth --bs 8
```

To apply SOM ():

```
python som.py
```

3.3 Data analysis

The Data Analysis Module provides a robust suite of tools for differential expression analysis and various downstream bioinformatics analyses. This module includes commonly used and advanced visualisation techniques such as:

- Correlation heatmaps
- Principal Component Analysis (PCA)
- Volcano plots integrated with Gene Set Enrichment Analysis (GSEA)
- Differential gene cluster heatmaps
- Venn diagrams in three styles: classic Venn, VennPie, and barplot
- Functional enrichment analysis, including:
 - Gene Ontology (GO) enrichment visualised in five styles: bubble, barplot, chord, cluster, and circle
 - KEGG pathway enrichment visualised in three formats: bubble, chord, and cluster
- Dumbbell-bar hybrid plots for comparative expression analysis

Additionally, this module provides flexible image export options suitable for

publication and presentation purposes.

Preparatory Work

Before running any analysis:

1. Set up the R environment as described in the "Environment" section.
2. Activate the R conda environment with:

```
conda activate r_env_name
```

(Replace 'r_env_name' with your environment name.)

3. Navigate to the analysis directory:

```
cd code/analysis
```

1: Correlation Heatmap

Script: 'cor_heatmap.r'

This script generates a correlation heatmap from the input expression matrix

Required Parameters:

- -i (input): Path to the input data file (tab-separated .txt file). The first row should contain column names.

Example Commands:

- Run with default settings:

```
Rscript cor_heatmap.r
```

- Run with a specified input file:

```
Rscript cor_heatmap.r -i cor_heatmap_test.txt
```




Figure 2. An example of a correlation heatmap.

2: Differential Gene Cluster Heatmap

Script: “df_cluster_heatmap.R”

This script generates a clustered heatmap of differentially expressed genes with optional row and column annotations.

Script Parameters:

- **-i (input):** Path to the input data file. The first row should contain column names, and the first column should contain row names. The file should be tab-delimited (t).
- **-a (type):** Specifies the type of figure to generate, including row and/or column annotations. Choose one of the following:
 - ⇒ **data_with_row_col:** Includes both row and column annotations.
 - ⇒ **data_with_col_annotation:** Includes only column annotations (e.g., Group, Age, Grade, Stage, Sex).
 - ⇒ **data_with_row_annotation:** Includes only row annotations (e.g., Path, Celltype).
 - ⇒ **'only_data':** indicates only output heatmap with no annotations.
- **-b (show_row_name):** Whether to display row names. Input ‘TRUE’ or

‘FALSE’.

- -c (show_col_name): Whether to display column names. Input ‘TRUE’ or ‘FALSE’.
- -d (clustering_dis_row): Method used for calculating row clustering distance. Choose from: ‘correlation’, ‘euclidean’, ‘maximum’, ‘manhattan’, ‘canberra’, ‘binary’, ‘minkowski’.
- -e (clustering_dis_col): Method used for calculating column clustering distance. Same options as above.
- -g (annotation_col_file): Path to the column annotation file. Required if ‘-a’ is set to ‘data_with_row_col’ or ‘data_with_col_annotation’.
- -f (annotation_row_file): Path to the row annotation file. Required if ‘-a’ is set to ‘data_with_row_col’ or ‘data_with_row_annotation’.
- -h (scal): Whether to scale the data. Choose from:
 - ⇒ row: Scale by row
 - ⇒ column: Scale by column
 - ⇒ none: No scaling
- -k (group): Whether to display grouping in the heatmap. Requires the column annotation file to include a ‘group’ column. Input ‘TRUE’ or ‘FALSE’.

Example Commands

Run the script with default parameters:

```
Rscript df_cluster_heatmap.R
```

Run the script with specified parameters:

```
Rscript df_cluster_heatmap.R -i df_gene_cluster_test.txt -a  
data_with_col_annotation -g df_gene_cluster_annotation_col.txt -c FALSE -b  
TRUE -d correlation -e correlation -h row -k TRUE
```




Figure 3. Example output of a differential gene cluster heatmap.

3: PCA

Script Parameters

- **-i (input):** the path of input data file. The first column of the data is the group information. The separator is '\t'.
- **-c (confidence_level):** confidence level.
- **-b (boundary):** decides whether to add boundary plots, which are the top and right sub-plots of the figure. Please input TRUE or FALSE.
- **-p (permanova):** decides whether to add PERMANOVA analysis. Please input TRUE or FALSE.
- **-m (method):** the method for permanova analysis. Please input one of 'manhattan', 'euclidean', 'canberra', 'clark', 'bray', 'kulczynski', 'jaccard', 'gower', 'altGower', 'morisita', 'horn', 'mountford', 'raup', 'binomial', 'chao',

'cao', 'mahalanobis', 'chisq', 'chord', 'hellinger', 'aitchison', and 'robust.aitchison'.

Example Commands

Run the script with default parameters:

```
Rscript pca.R
```

Run the script with specified parameters:

```
Rscript pca.R -i pca_test.txt -c 0.92 -b FALSE -p FALSE
```




Figure 4. An example of PCA.

4: VENN

Script Parameters

- -i (input): the path of input data file. The first row is the column names, and each column is the information for each group. The separator is '\t'.
- -t (type): the type of plot. Please input one of Venn, Vennpie and Barplot.

Example Command

Run the script with default parameters:

```
Rscript venn.R
```

Run the script with specified parameters:

```
Rscript venn.R -i venn_test.txt -t Vennpie
```




Figure 5. An example of VENN.

5: Volcano with GSEA

Script Parameters

- -i (input): the path of input data file. Keep logFC, padj, gene column names of dataset consistent. The separator is '\t'.
- -g (gmt): the path of gmt data file. If you want to conduct GSEA analysis, then specialize it to the path gmt data file, otherwise set it to 'none'.
- -a (fc_thr): the threshold of logFC.
- -b (padj_thr): the threshold of padj (adjusted p-value).
- -c (top): the number of higher-order genes to be emphasized. The emphasized genes are in the dashed box.
- -d (top_fc_thr): the threshold of logFC for the higher-order gene. It should be more strict than 'fc_thr'.

- -e (top_padj_thr): the threshold of padj for the higher-order gene. It should be more strict than 'padj_thr'.
- -f (gene_sig): the list of marked genes, the separator must be comma ','.

Example Commands

Run the script with default parameters:

```
Rscript volcano_gsea_padj.R
```

Run the script with specified parameters:

```
Rscript volcano_gsea_padj.R -i volcano_test.txt -g none -a 1 -b 0.05 -c 30 -d 2 -e 0.01 -f KRAS,FOSL1
```




Figure 6. An example of volcano with GSEA.

6: GO Enrichment

Script Parameters

- -i (input): the path of input data file. The first column is gene symbol, and the column name is gene. The data file only needs the gene symbol if 'type' is bubble or bar, otherwise needs to make sure that the data file has numeric data with the column name logFC to sort by the size of the value.

- -a (type): the type of figure. Please input one of bubble, bar, chord, cluster and circle.
- -b (organism): Please input one of Homo_sapiens, Mus_musculus, Bovine, Homo_sapiens and Drosophila_melanogaster.
- -c (pvalue): pvalue threshold for GO enrichment analysis.
- -d (qvalue): qvalue threshold for GO enrichment analysis.
- -e (termNum): the number of terms for each ontology to be displayed.

Example Commands

Run the script with default parameters:

```
Rscript go_enrichment.R
```

Run the script with specified parameters:

```
Rscript go_enrichment.R -i go_enrichment_test.txt -a bubble -b Mus_musculus -c 0.01 -e 10
```




Figure 7. An example of GO Enrichment.

7: KEGG Enrichment

Script Parameters

- -i (input): the path of input data file. The first column is gene symbol, and the column name is gene. The data file only needs the gene symbol if

'type' is bubble, otherwise needs to make sure that the data file has numeric data with the column name logFC to sort by the size of the value.

- -a (type): the type of figure. Please input one of bubble, chord, and cluster.
- -b (organism): Please input one of hsa, mmu, bos, and dme.
- -c (pvalue): pvalue threshold for KEGG enrichment analysis.
- -d (qvalue): qvalue threshold for KEGG enrichment analysis.
- -e (termNum): the number of terms for each ontology to be displayed.

Example Commands

Run the script with default parameters:

```
Rscript kegg_enrichment.R
```

Run the script with specified parameters:

```
Rscript kegg_enrichment.R -i kegg_enrichment_test.txt -a bubble -b dme -d 0.01 -e 10
```




Figure 8. An example of KEGG Enrichment.

8: Dumbbell_Bar

Script Parameters

- -i (input): the path of input data file. This file is used to draw dumbbell diagrams.
- -c (data_bar): the file path of bar plot. This file is used to draw barplot

diagrams. Keep the contents and name of the first column in both two files similar and same respectively.

- -a (y_label): the y label of dumbbell_bar plot.
- -b (mark_fams): the terms that will be marked, which are factors in the first column of data files, the separator must be comma ','.

Example Commands

Run the script with default parameters:

```
Rscript dumbbell_bar.R
```

Run the script with specified parameters:

```
Rscript dumbbell_bar.R -i dumbbell_test.txt -c dumbbell_barplot_test.txt -a number  
-b Notothenioid
```




Figure 9. An example of Dumbbell_Bar.

9: DESeq2 for Read Count

Script Parameters

- -i (expression_file): the path of expression file. The first column is row names, the first row is column names.
- -k (info_file): the path of group information file. You need to make sure that the order of the samples in the expression file corresponds to the group info file order here. Keep the row names in the group file the same as the column names in the expression file: Control_1, Control_2, Treatment_1, Treatment_2. The Group file must contain a Group column, and the value of the group column must be 'Control' or 'Treatment'.
- -c (fc_thr): the log2FC threshold for differential expression analysis.
- -d (padj_thr): the padj threshold for differential expression analysis.

Example Commands

Run the script with default parameters:

```
Rscript DESeq2_read_count.R
```

Run the script with specified parameters:

```
Rscript DESeq2_read_count.R -i expression.txt -k info.txt -c 2 -d 0.05
```

10: limma for FPKM

Script Parameters

- -i (expression_file): the path of expression file. The first column is row names, the first row is column names.
- -k (info_file): the path of group information file. You need to make sure that the order of the samples in the expression file corresponds to the group info file order here. Keep the row names in the group file the same as the column names in the expression file: Control_1, Control_2, Treatment_1, Treatment_2. The Group file must contain a Group column, and the value of the group column must be 'Control' or 'Treatment'.

- -c (fc_thr): the log2FC threshold for differential expression analysis.
- -d (padj_thr): the padj threshold for differential expression analysis.

Example Commands

Run the script with default parameters:

```
Rscript limma_fpkm_df.R
```

Run the script with specified parameters:

```
Rscript limma_fpkm_df.R -i expression.txt -k info.txt -c 2 -d 0.05
```